



Scan 4 all  
introduces



INSIDE / **My partners**

# A CONNECTED ENVIRONMENT

IRIS BY STARCK® allows you to collaborate with different types of partners to facilitate the management of your orders and communications.



## Depending on your needs, you can :

- › invite your own contacts,
- › connect a partner application,
- › automatically sync your data with certain compatible services.

## How do I connect with my partners?

### Via IRIS CONNECT (invitez vos propres contacts)

Contacts refer to the partners you wish to add/invite via the IRIS CONNECT directory.

Can't find your contact?

Send an invitation by email


contact@mail.fr

Invite

From IRIS CONNECT, use the "search" field to quickly find a collaborator. If they are not part of the IRIS directory, enter their email address to invite them. They will receive a registration link associated with a unique security code to create their web account at [www.irisbystarck.com](http://www.irisbystarck.com).

Once created, these contacts can receive your orders and interact with you directly via messages in IRIS CONNECT. They can also access the dedicated web interface, allowing them to view, track, and process requests sent to [www.irisbystarck.com](http://www.irisbystarck.com).


### Via pairing (connect with our partner software)



If you have an account with one of the partners integrated with IRIS BY STARCK® (e.g., SMILERS, CIRCLE, etc.): **you can send your order directly to your account** with that partner.

**When creating an order and selecting a partner**, click on that partner's card and follow the prompts to complete the pairing.

## Via data synchronization (LaGalaxy by BIOTECH DENTAL)



All your order information is automatically synchronized between IRIS BY STARCK® and the partner service—such as LaGalaxy—ensuring continuous data updates. More information at


[www.biotechgalaxy.com](http://www.biotechgalaxy.com)

LaGalaxy  
By Biotech Dental




## How to delete a partner?

### From the IRIS CONNECT directory




You can delete a contact directly from **IRIS CONNECT > My Contacts** by clicking the “trash can” icon on their contact card.

### For paired third-party partner software



Go to **My Account > Manage Partner Accounts** and click “Forget” to revoke the pairing of a partner service.

### For synchronized software (LaGalaxy)



Synchronization can be disabled at any time from **My Account** by accessing the associated partner service.

Unlinking stops data synchronization and revokes the permissions you granted to the service.

# GO FURTHER WITH IRIS BY STARCK®.

Explore all the resources available to you.



Inside / [Quick start](#)



Inside / [My settings](#)



Inside / [DIRECT SCAN](#)



Inside / [My account](#)



Inside / [System informations](#)



Inside / [Tools related to scanning](#)



Inside / [IRIS Connect](#)



Inside / [SAV](#)



Inside / [My orders](#)



Inside / [Real time scan](#)



Inside / [My notifications](#)



Inside / [Computer specifications](#)



Inside / [My partners](#)




Inside / [Guided workflow](#)



Inside / [My patients](#)



Inside / [Open workflow](#)



Access all the  
documentation.

## Manufacturer

**Scan<sup>4</sup>all**

Excellence, for everyone.

137 rue Garance 31670 Labège - FRANCE

N° SIRET 507 484 863 00059

[www.linkedin.com/company/scan4all/](https://www.linkedin.com/company/scan4all/)

IRIS BY STARCK®

Class I medical device (MD) intended for digital dentistry. CE.

For healthcare professionals only.

© Scan4all - Version 1.0 - 2026